pjama.eu

GB

INTRODUCTION

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PJAMA SENSOR

INSTRUCTION

For best results: The child should be motivated and want to get dry. The parent needs to help the child wake up, especially at the beginning of treatment.

NOTE: The sensor does not need to be charged. Its equipped with Bluetooth® and the sensor communicates directly with the **Pjama Speaker** or **DryGuardians SmartWatch**.

1. Turn On your Pjama Speaker.
2. Test the Sensor with the Audio Device. Short-circuit the Sensor buttons by using a metallic object. This will trigger the alarm. The sensor range is 10 meters (30 feet).
3. Attach the sensor to the buttons located on the inside of the Pjama Treatment Pants/ Shorts or on the outside of the Underwear.
4. The sensor is now ready for use.
5. Disconnect the sensor from the Treatment Garment **immediately** after the alarm has sounded. It prevents the alarm from sounding again and saves the battery.

The Sensor can be cleaned with soap and water or disinfected with alcohol-based disinfectants.

WARRANTY

The Pjama Sensor has a one- year warranty from the buyer's purchase date of the alarm. You must be the original buyer and show evidence of your purchase. This warranty does not apply if malfunction is caused by an accident or similar incident or injury, water damage, negligence, abnormal use or failure to comply with the instructions in this manual.

CAUTION

- Never place the sensor directly against the skin.
- Pjama Sensor is intended for use as part of an enuresis alarm.
- Make sure that no part of the Sensor is put in the mouth for any reason. This could cause choking or other harm.
- The Pjama Sensor is a medical device that should be treated with care.
- Do not hold an alarming Speaker /Audio Device within 10 cm (4 inches) of the ear. This could cause damage to the ear.
- Disconnect the sensor from the Treatment Garment immediately after the alarm has sounded. It saves the battery time on the sensor.
- Avoid storing the Sensor on metal surfaces to prevent interference or damage.

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